

Transform your endodontic practice with
Statera vi, where precision meets innovation


Cordless Endo Motor with Integrated Apex Locator

Statera vi empowers you to achieve unparalleled results with Forumtec's patented **Balanced Force adaptive technique**, elevating to a new level of excellence



STATERA vi

Rotating 360° Contra Angle
for improved control and
ease of handling



High Resolution
Color Screen

High speed settings
ranging from 50 to
3200 RPM and torque
levels reaching up to
6.0 Ncm.



STATERA  **vi**

Precise measurements with Forumtec's integrated apex locator patented technology



Magnetic cradle with embedded power bank for automated wireless charging

Reciprocation motion with a 5° interval and adjustment range from 5° to 400°



STATERA  **vi**

Advanced features that empower you to achieve unparalleled results with confidence



Efficiency and Enhanced Safety

Utilizing Forumtec's patented technology, Statera vi provides unparalleled versatility, enabling you to utilize a wide range of files in Balanced Force reciprocation mode for efficient and safe canal preparation



Increased durability

Brushless motor for superior performance ensuring a longer lifespan, and reduced maintenance



Ergonomic design

The triangular shape of the hand piece is designed for ergonomic comfort, ensuring stable placement on the tray.

